

## Antibodies Product Datasheet

<b>Product Name:</b>	NDRG1 (Phospho-Thr346) Antibody
<b>Catalog Number:</b>	A9595
<b>Reactivity:</b>	Human, Mouse, Rat
<b>Applications:</b>	WB
<b>Dilutions:</b>	WB 1:1000-2000
<b>Modification:</b>	Phospho Specific
<b>Source:</b>	Polyclonal Rabbit
<b>Purification:</b>	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
<b>Immunogen:</b>	Synthesized phospho peptide around human NDRG1 (Thr346)
<b>Concentration:</b>	1 mg/ml
<b>Gene Name:</b>	NDRG1 CAP43 DRG1 RTP
<b>Protein Name:</b>	NDRG1 (Thr346)
<b>Other Name:</b>	Protein NDRG1 (Differentiation-related gene 1 protein) (DRG-1) (N-myc downstream-regulated gene 1 protein) (Nickel-specific induction protein Cap43) (Reducing agents and tunicamycin-responsive protein) (RTP) (Rit42)
<b>Human Swiss Prot No.:</b>	Q92597
<b>Formulation:</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Specificity:</b>	This antibody detects endogenous levels of Human Mouse Rat NDRG1 (phospho-Thr346)
<b>Storage and Stability:</b>	-20°C/1 year