

Antibodies Product Datasheet

Product Name:	SIRT6 mouse mAb
Catalog Number:	N0839
Reactivity:	
Applications:	WB
Dilutions:	wb 1:500
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Purified recombinant human SIRT6 protein fragments expressed in E.coli.
Concentration:	1 mg/ml
Gene Name:	sirt6
Protein Name:	
Other Name:	2810449N18Rik;AI043036;Mono ADP ribosyltransferase sirtuin 6;NAD-dependent deacetylase sirtuin-6;NAD-dependent protein deacetylase sirtuin-6;Regulatory protein SIR2 homolog;SIR2 like 6;SIR2 like protein 6;Sir2 related protein type 6;SIR2-like protein 6;SI
Human Swiss Prot No.:	Q8N6T7
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Transfected Only.
Storage and Stability:	-20°C/1 year