

Antibodies Product Datasheet

Product Name:	C/EBP α Antibody
Catalog Number:	B0056
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity- chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human C/EBP-alpha. AA range:6-55
Concentration:	1 mg/ml
Gene Name:	СЕВРА
Protein Name:	CCAAT/enhancer-binding protein alpha
Other Name:	CEBPA; CCAAT/enhancer-binding protein alpha; C/EBP alpha
Human Swiss Prot No.:	P49715
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	C/EBP α Polyclonal Antibody detects endogenous levels of C/EBP α protein.
Storage and Stability:	-20°C/1 year