

Assay Biotechnology Company

47787 Fremont Blvd. Fremont, CA 94538 Phone: (408) 747 0185 Fax: (408 747 0145 www.assaybiotechnlogy.com

Antibodies Product Datasheet

Product Name: EpoR Antibody

Catalog Number: B0917

Reactivity: Human, Mouse, Rat, Monkey

Applications: WB, IF/ICC, ELISA

Dilutions: Western Blot: 1/500 - 1/2000. Immunofluorescence: 1/200 - 1/1000.

ELISA: 1/40000. Not yet tested in other applications.

Modification: Unmodified/Total

Source: Polyclonal Rabbit

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-

chromatography using epitope-specific immunogen.

Immunogen: The antiserum was produced against synthesized peptide derived from

human Epo-R. AA range:341-390

Concentration: 1 mg/ml

Gene Name: EPOR

Protein Name: Erythropoietin receptor

Other Name: EPOR; Erythropoietin receptor; EPO-R

Human Swiss Prot No.: P19235

Formulation: Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium

azide.

Specificity: EpoR Polyclonal Antibody detects endogenous levels of EpoR protein.

Storage and Stability: -20°C/1 year