

Antibodies Product Datasheet

Product Name:	Anti-Human Cytokeratin 8 (ABT537) Mouse Monoclonal Antibody
Catalog Number:	V0494
Clone ID:	ABT537
Reactivity:	Human, Mouse, Rat
Applications:	IHC, WB
Application Dilutions	IHC 1:100-200, WB 1:500-1000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
UniProt:	KRT8
Gene Name:	P05787
Other Name:	CARD2, CK 8, CK-8, CK8, CYK8, CYKER, Cytokeratin endo A, Cytokeratin-8, DreK8, EndoA, K2C8, K2C8_HUMAN, K8, Keratin 8, Keratin type II cytoskeletal 8, Keratin, type II cytoskeletal 8, Keratin-8, KO, Krt 2.8, KRT8, MGC118110, MGC174782, MGC53564, MGC85764, sb:cb186, Type-II keratin Kb8
Isotype	IgG2b, Kappa
Localization:	Liver/ Tonsil
Positive Tissue Control:	Cytoplasmic, Membranous
Immunogen	Synthetic peptide
Purification	Immunogen affinity purified
Formulation:	PBS, pH7.4, 50% glycerol, 0.03%Proclin300
Specificity:	The antibody can specifically recognize human CK8 protein, and shows no cross reaction with CK1, 4, 5, 6, 7, 10, 15, 16, 17, 18, 20.
Storage and Stability:	-20°C/1 year