

## Antibodies Product Datasheet

<b>Product Name:</b>	p27 Antibody
<b>Catalog Number:</b>	B0965
<b>Reactivity:</b>	Human, Rat, Mouse
<b>Applications:</b>	WB, IHC-p, IF(paraffin section), ELISA
<b>Dilutions:</b>	WB 1:500-2000, ELISA 1:10000-20000 IHC 1:50-300
<b>Modification:</b>	Unmodified/Total
<b>Source:</b>	Polyclonal Rabbit
<b>Purification:</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Immunogen:</b>	The antiserum was produced against synthesized peptide derived from human p27 Kip1. AA range:149-198
<b>Concentration:</b>	1 mg/ml
<b>Gene Name:</b>	CDKN1B
<b>Protein Name:</b>	Cyclin-dependent kinase inhibitor 1B
<b>Other Name:</b>	CDKN1B; KIP1; Cyclin-dependent kinase inhibitor 1B; Cyclin-dependent kinase inhibitor p27; p27Kip1
<b>Human Swiss Prot No.:</b>	P46527
<b>Formulation:</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Specificity:</b>	p27 Polyclonal Antibody detects endogenous levels of p27 protein.
<b>Storage and Stability:</b>	-20°C/1 year