
Antibodies Product Datasheet

Product Name:	TEL Antibody
Catalog Number:	B8432
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human Tel. AA range:223-272
Concentration:	1 mg/ml
Gene Name:	ETV6
Protein Name:	Transcription factor ETV6
Other Name:	ETV6; TEL; TEL1; Transcription factor ETV6; ETS translocation variant 6; ETS-related protein Tel1; Tel
Human Swiss Prot No.:	P41212
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	TEL Polyclonal Antibody detects endogenous levels of TEL protein.
Storage and Stability:	-20°C/1 year