

Antibodies Product Datasheet

Product Name:	Cyclin H mouse mAb
Catalog Number:	N0855
Reactivity:	Human
Applications:	WB,IP
Dilutions:	wb 1:1000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Purified recombinant human Cyclin H protein fragments expressed in E.coli.
Concentration:	1 mg/ml
Gene Name:	ccnh
Protein Name:	
Other Name:	6330408H09Rik;AI661354;AV102684;AW538719;CAK;CAK complex subunit;ccnh;CCNH_HUMAN;CDK activating kinase;CDK activating kinase complex subunit;Cyclin dependent kinase activating kinase;cyclin dependent kinase activating kinase complex subunit;Cyclin H; Cycl
Human Swiss Prot No.:	P51946
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Cyclin H and does not cross-react with related proteins.
Storage and Stability:	-20°C/1 year