
Antibodies Product Datasheet

Product Name:	IFT20 Antibody
Catalog Number:	C44217
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human IFT20
Concentration:	1 mg/ml
Gene Name:	IFT20
Protein Name:	Intraflagellar transport protein 20 homolog (hIFT20)
Other Name:	
Human Swiss Prot No.:	Q8IY31
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human IFT20
Storage and Stability:	-20°C/1 year