

Antibodies Product Datasheet

Product Name:	IGFBP1 Antibody
Catalog Number:	C30203
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. IHC-p: 1:100-1:300. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from the Internal region of human IGFBP1. AA range:91-140
Concentration:	1 mg/ml
Gene Name:	IGFBP1
Protein Name:	Insulin-like growth factor-binding protein 1
Other Name:	IGFBP1; IBP1; Insulin-like growth factor-binding protein 1; IBP-1; IGF-binding protein 1; IGFBP-1; Placental protein 12; PP12
Human Swiss Prot No.:	P08833
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	IGFBP1 Polyclonal Antibody detects endogenous levels of IGFBP1 protein.
Storage and Stability:	-20°C/1 year