

Antibodies Product Datasheet

Product Name:	CKR-6 Antibody
Catalog Number:	C30470
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from the N-terminal region of human CCR6. AA range:1-50
Concentration:	1 mg/ml
Gene Name:	CCR6
Protein Name:	C-C chemokine receptor type 6
Other Name:	CCR6; CKRL3; CMKBR6; GPR29; STRL22; C-C chemokine receptor type 6; C-C CKR-6; CC-CKR-6; CCR-6; Chemokine receptor-like 3; CKR-L3; DRY6; G-protein coupled receptor 29; GPR-CY4; GPRCY4; LARC receptor; CD196
Human Swiss Prot No.:	P51684
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CKR-6 Polyclonal Antibody detects endogenous levels of CKR-6 protein.
Storage and Stability:	-20°C/1 year