

Antibodies Product Datasheet

Product Name:	SUMF1 Antibody
Catalog Number:	C50651
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human SUMF1
Concentration:	1 mg/ml
Gene Name:	SUMF1 FGE PSEC0152 UNQ3037/PRO9852
Protein Name:	Sulfatase-modifying factor 1 (EC 1.8.99.-) (C-alpha-formylglycine-generating enzyme 1)
Other Name:	
Human Swiss Prot No.:	Q8NBK3
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of SUMF1 at Human, Mouse
Storage and Stability:	Store at -20° C/1 year