

Antibodies Product Datasheet

Product Name:	C1S (light chain, Cleaved-Ile438) Antibody
Catalog Number:	L0406
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:1000-2000 ELISA 1:5000-20000
Modification:	Cleaved Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human C1S (light chain, Cleaved-Ile438)
Concentration:	1 mg/ml
Gene Name:	C1S
Protein Name:	C1S (light chain, Cleaved-Ile438)
Other Name:	Complement C1s subcomponent (EC 3.4.21.42) (C1 esterase) (Complement component 1 subcomponent s) [Cleaved into: Complement C1s subcomponent heavy chain; Complement C1s subcomponent light chain]
Human Swiss Prot No.:	P09871
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human C1S (light chain, Cleaved-Ile438, protein was cleaved amino acid sequence between 437-438)
Storage and Stability:	-20°C/1 year