

Antibodies Product Datasheet

Product Name:	JAK3 Mouse mAb(6B5)
Catalog Number:	N1485
Reactivity:	Human, Mouse,Rat
Applications:	IHC
Dilutions:	IHC 1:50-200
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human JAK3
Concentration:	1 mg/ml
Gene Name:	JAK3
Protein Name:	JAK3
Other Name:	Tyrosine-protein kinase JAK3 (EC 2.7.10.2) (Janus kinase 3) (JAK-3) (Leukocyte janus kinase) (L-JAK)
Human Swiss Prot No.:	P52333
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.58% sodium azide.
Specificity:	This antibody detects endogenous levels of JAK3 at Human, Mouse,Rat
Storage and Stability:	-20°C/1 year