

## Antibodies Product Datasheet

<b>Product Name:</b>	Cdk4 Monoclonal Antibody
<b>Catalog Number:</b>	N0673
<b>Reactivity:</b>	Human, Mouse, Rat, Dog, Pig
<b>Applications:</b>	WB
<b>Dilutions:</b>	Western Blot: 1/1000 - 1/2000. Not yet tested in other applications.
<b>Modification:</b>	Unmodified/ Total
<b>Source:</b>	Monoclonal Mouse
<b>Purification:</b>	Affinity purification
<b>Immunogen:</b>	Purified recombinant human Cdk4 protein fragments expressed in E.coli.
<b>Concentration:</b>	1 mg/ml
<b>Gene Name:</b>	CDK4
<b>Protein Name:</b>	Cell division protein kinase 4
<b>Other Name:</b>	CDK4; Cyclin-dependent kinase 4; Cell division protein kinase 4; PSK-J3
<b>Human Swiss Prot No.:</b>	P11802
<b>Formulation:</b>	Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50% glycerol.
<b>Specificity:</b>	Cdk4 Monoclonal Antibody detects endogenous levels of Cdk4 protein.
<b>Storage and Stability:</b>	-20°C/1 year