

Antibodies Product Datasheet

Product Name:	TREM-1 Antibody
Catalog Number:	C22137
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. IHC-p: 1:100-300 ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from TREM-1 . at AA range: 40-120
Concentration:	1 mg/ml
Gene Name:	TREM1
Protein Name:	Triggering receptor expressed on myeloid cells 1
Other Name:	TREM1; Triggering receptor expressed on myeloid cells 1; TREM-1; Triggering receptor expressed on monocytes 1; CD354
Human Swiss Prot No.:	Q9NP99
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	TREM-1 Polyclonal Antibody detects endogenous levels of TREM-1 protein.
Storage and Stability:	-20°C/1 year