

Antibodies Product Datasheet

Product Name:	ECEL1 Antibody
Catalog Number:	C43571
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human ECEL1
Concentration:	1 mg/ml
Gene Name:	ECEL1 XCE UNQ2431/PRO4991
Protein Name:	ECEL1
Other Name:	Endothelin-converting enzyme-like 1 (EC 3.4.24.-) (Xce protein)
Human Swiss Prot No.:	O95672
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.360% sodium azide.
Specificity:	This antibody detects endogenous levels of ECEL1 at Human/Mouse/Rat
Storage and Stability:	-20°C/1 year