

Antibodies Product Datasheet

Product Name:	CDK5 Antibody
Catalog Number:	R12-2605
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC, ELISA
Dilutions:	CDK5 Antibody
Modification:	Unmodified/Total
Source:	Rabbit Polyclonal
Purification:	The antibody was purified from rabbit antiserum by affinity-chromatography using immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human CDK5.
Concentration:	1 mg/ml
Gene Name:	CDK5
Other Name:	CDK5; Cell division protein kinase 5; CRK6; kinase Cdk5; PSSALRE; Tau protein kinase II catalytic subunit; TAU protein kinase II catalytic subunit; TPKII catalytic subunit
Human Swiss Prot No.:	Q00535
Formulation:	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Specificity:	CDK5 Antibody antibody detects endogenous levels of CDK5.
Storage and Stability:	Stable at -20°C for at least 1 year.