

Antibodies Product Datasheet

Product Name:	TCF-3 Antibody
Catalog Number:	C10385
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human TCF3. AA range:571-620
Concentration:	1 mg/ml
Gene Name:	TCF3
Protein Name:	Transcription factor E2-alpha
Other Name:	TCF3; BHLHB21; E2A; ITF1; Transcription factor E2-alpha; Class B basic helix-loop-helix protein 21; bHLHb21; Immunoglobulin enhancer-binding factor E12/E47; Immunoglobulin transcription factor 1; Kappa-E2-binding factor; Transcription facto
Human Swiss Prot No.:	P15923
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	TCF-3 Polyclonal Antibody detects endogenous levels of TCF-3 protein.
Storage and Stability:	-20°C/1 year