

Antibodies Product Datasheet

Product Name:	MAGT1 Antibody
Catalog Number:	C48130
Reactivity:	Human, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human MAGT1
Concentration:	1 mg/ml
Gene Name:	MAGT1 IAG2 PSEC0084 UNQ628/PRO1244
Protein Name:	Magnesium transporter protein 1 (MagT1) (Implantation-associated protein) (IAP)
Other Name:	
Human Swiss Prot No.:	Q9H0U3
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of MAGT1 at Human,Rat
Storage and Stability:	Store at -20° C/1 year