

Antibodies Product Datasheet

Product Name:	Caspase-8 Monoclonal Antibody
Catalog Number:	N0092
Reactivity:	Human, Mouse, Rat, Monkey
Applications:	WB, IHC-p, IF (paraffin section), FCM, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/200 - 1/1000. Flow cytometry: 1/200 - 1/400. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of human Caspase-8 expressed in E. Coli.
Concentration:	
Gene Name:	CASP8
Protein Name:	Caspase8
Other Name:	CASP8; MCH5; Caspase-8; CASP-8; Apoptotic cysteine protease; Apoptotic protease Mch-5; CAP4; FADD-homologous ICE/ced-3-like protease; FADD-like ICE; FLICE; ICE-like apoptotic protease 5; MORT1-associated ced-3 homolog; MACH
Human Swiss Prot No.:	Q14790
Formulation:	Ascitic fluid containing 0.03% sodium azide, 0.5% BSA, 50% glycerol.
Specificity:	Caspase-8 Monoclonal Antibody detects endogenous levels of Caspase-8 protein.
Storage and Stability:	-20°C/1 year