

## Antibodies Product Datasheet

<b>Product Name:</b>	Histone H1 (Phospho-Thr17) Antibody
<b>Catalog Number:</b>	A9710
<b>Reactivity:</b>	Human, Mouse, Rat
<b>Applications:</b>	WB
<b>Dilutions:</b>	WB: 1:1000-2000
<b>Modification:</b>	Phospho Specific
<b>Source:</b>	Polyclonal Rabbit
<b>Purification:</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using specific immunogen.
<b>Immunogen:</b>	Synthetic Peptide of Histone H1 (Phospho Thr17)
<b>Concentration:</b>	
<b>Gene Name:</b>	HIST1H1A/HIST1H1C/HIST1H1D/HIST1H1E
<b>Protein Name:</b>	Histone H1.1/Histone H1.2/Histone H1.3/Histone H1.4
<b>Other Name:</b>	HIST1H1A; H1F1; Histone H1.1; Histone H1a; HIST1H1C; H1F2; Histone H1.2; Histone H1c; Histone H1d; Histone H1s-1; HIST1H1D; H1F3; Histone H1.3; Histone H1c; Histone H1s-2; HIST1H1E; H1F4; Histone H1.4; Histone H1b; Histone H1s-4; YM3309
<b>Human Swiss Prot No.:</b>	Q02539/P16403/P16402/P10412
<b>Formulation:</b>	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
<b>Specificity:</b>	The antibody detects endogenous Histone H1 (Phospho Thr17) protein.
<b>Storage and Stability:</b>	-20°C/1 year