

Antibodies Product Datasheet

Product Name:	ASAH1 Antibody
Catalog Number:	C48293
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human ASAH1
Concentration:	1 mg/ml
Gene Name:	ASAH1 ASAH HSD-33 HSD33
Protein Name:	Acid ceramidase (AC) (ACDase) (Acid CDase) (EC 3.5.1.23) (Acylsphingosine deacylase) (N-acylsphingosine amidohydrolase) (Putative 32 kDa heart protein) (PHP32) [Cleaved into: Acid ceramidase subunit alphaAcid ceramidase subunit beta]
Other Name:	
Human Swiss Prot No.:	Q13510
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of ASAH1 at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year