

Antibodies Product Datasheet

Product Name:	MDM2 (Phospho-Ser186/S188) Antibody
Catalog Number:	A1154
Reactivity:	Human, Rat, Mouse
Applications:	IF/ICC, WB, ELISA
Dilutions:	WB 1:500-2000 Immunofluorescence: 1/200 - 1/1000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human MDM2 around the phosphorylation site of Ser186 and Ser188. AA range:151-200
Concentration:	1 mg/ml
Gene Name:	MDM2
Protein Name:	E3 ubiquitin-protein ligase Mdm2
Other Name:	MDM2; E3 ubiquitin-protein ligase Mdm2; Double minute 2 protein; Hdm2; Oncoprotein Mdm2; p53-binding protein Mdm2
Human Swiss Prot No.:	Q00987
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Phospho-MDM2 (S186/S188) Polyclonal Antibody detects endogenous levels of MDM2 protein only when phosphorylated at S186/S188.
Storage and Stability:	-20°C/1 year