
Antibodies Product Datasheet

Product Name:	Tel Antibody
Catalog Number:	R12-3611
Reactivity:	Human
Applications:	WB, ELISA
Dilutions:	Tel Antibody
Modification:	Unmodified/Total
Source:	Rabbit Polyclonal
Purification:	The antibody was purified from rabbit antiserum by affinity-chromatography using immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human Tel.
Concentration:	1 mg/ml
Gene Name:	ETV6
Other Name:	ETS translocation variant 6; ets variant gene 6 (TEL oncogene); ETS-related protein Tel1; ETV6; TEL1; Transcription factor ETV6
Human Swiss Prot No.:	P41212
Formulation:	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Specificity:	Tel Antibody antibody detects endogenous levels of Tel.
Storage and Stability:	Stable at -20°C for at least 1 year.