

Antibodies Product Datasheet

Product Name:	cPLA2- ϵ Antibody
Catalog Number:	C15307
Reactivity:	Human, Rat, Mouse
Applications:	IHC-p, IF(paraffin section), ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human PLA2G4E. AA range:401-450
Concentration:	1 mg/ml
Gene Name:	PLA2G4E
Protein Name:	Cytosolic phospholipase A2 epsilon
Other Name:	PLA2G4E; Cytosolic phospholipase A2 epsilon; cPLA2-epsilon; Phospholipase A2 group IVE
Human Swiss Prot No.:	Q3MJ16
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	cPLA2- ϵ Polyclonal Antibody detects endogenous levels of cPLA2- ϵ protein.
Storage and Stability:	-20°C/1 year