

Antibodies Product Datasheet

Product Name:	TSEN54 Antibody
Catalog Number:	C19290
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human TSEN54. AA range:261-310
Concentration:	1 mg/ml
Gene Name:	TSEN54
Protein Name:	tRNA-splicing endonuclease subunit Sen54
Other Name:	TSEN54; SEN54; tRNA-splicing endonuclease subunit Sen54; SEN54 homolog; HsSEN54; tRNA-intron endonuclease Sen54
Human Swiss Prot No.:	Q7Z6J9
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	TSEN54 Polyclonal Antibody detects endogenous levels of TSEN54 protein.
Storage and Stability:	-20°C/1 year