
Antibodies Product Datasheet

Product Name:	CA IX Monoclonal Antibody
Catalog Number:	N0085
Reactivity:	Human
Applications:	WB,IHC-p,IF(paraffin section),FCM,ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/200 - 1/1000. Flow cytometry: 1/200 - 1/400. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of human CA IX expressed in E. Coli.
Concentration:	
Gene Name:	CA9
Protein Name:	Carbonic anhydrase 9
Other Name:	CA9; G250; MN; Carbonic anhydrase 9; Carbonate dehydratase IX; Carbonic anhydrase IX; CA-IX; CAIX; Membrane antigen MN; P54/58N; Renal cell carcinoma-associated antigen G250; RCC-associated antigen G250; pMW1
Human Swiss Prot No.:	Q16790
Formulation:	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Specificity:	CA IX Monoclonal Antibody detects endogenous levels of CA IX protein.
Storage and Stability:	-20°C/1 year