

Antibodies Product Datasheet

Product Name:	S-100A10 Monoclonal Antibody
Catalog Number:	N0572
Reactivity:	Human
Applications:	WB,IHC-p,IF/ICCELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/200 - 1/1000. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of human S-100A10 expressed in E. Coli.
Concentration:	
Gene Name:	S100A10
Protein Name:	Protein S100-A10
Other Name:	S100A10; ANX2LG; CAL1L; CLP11; Protein S100-A10; Calpactin I light chain; Calpactin-1 light chain; Cellular ligand of annexin II; S100 calcium-binding protein A10; p10 protein; p11
Human Swiss Prot No.:	P60903
Formulation:	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Specificity:	S-100A10 Monoclonal Antibody detects endogenous levels of S-100A10 protein.
Storage and Stability:	-20°C/1 year