
Antibodies Product Datasheet

Product Name:	LAT Antibody
Catalog Number:	R12-2226
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	LAT Antibody
Modification:	Unmodified/Total
Source:	Rabbit Polyclonal
Purification:	The antibody was purified from rabbit antiserum by affinity-chromatography using immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human LAT.
Concentration:	1 mg/ml
Gene Name:	LAT
Other Name:	Linker for activation of T-cells family member 1; 36 kDa phospho-tyrosine adapter protein
Human Swiss Prot No.:	O43561
Formulation:	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Specificity:	LAT antibody detects endogenous levels of LAT.
Storage and Stability:	Stable at -20°C for at least 1 year.