

Antibodies Product Datasheet

Product Name:	Flotillin-2 Antibody
Catalog Number:	R12-2143
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human Flotillin-2. AA range:95-144
Concentration:	1 mg/ml
Gene Name:	FLOT2
Protein Name:	Flotillin-2
Other Name:	FLOT2; ESA1; M17S1; Flotillin-2; Epidermal surface antigen; ESA; Membrane component chromosome 17 surface marker 1
Human Swiss Prot No.:	Q14254
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Flotillin-2 Polyclonal Antibody detects endogenous levels of Flotillin-2 protein.
Storage and Stability:	-20°C/1 year