

Antibodies Product Datasheet

Product Name:	Caspase-1 Antibody
Catalog Number:	B1025
Reactivity:	Human, Mouse, Rat
Applications:	IF/ICC, WB, IHC-p, ELISA
Dilutions:	IF: 1:50-200 WB 1:500-2000, IHC 1:50-300 IHC 1:50-300
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human Caspase 1. AA range:342-391
Concentration:	1 mg/ml
Gene Name:	CASP1
Protein Name:	Caspase1
Other Name:	CASP1; IL1BC; IL1BCE; Caspase-1; CASP-1; Interleukin-1 beta convertase; IL-1BC; Interleukin-1 beta-converting enzyme; ICE; IL-1 beta-converting enzyme; p45
Human Swiss Prot No.:	P29466
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Caspase-1 Polyclonal Antibody detects endogenous levels of Caspase-1 protein.
Storage and Stability:	-20°C/1 year