
Antibodies Product Datasheet

Product Name:	FES Antibody
Catalog Number:	R12-2768
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	FES Antibody
Modification:	Unmodified/Total
Source:	Rabbit Polyclonal
Purification:	The antibody was purified from rabbit antiserum by affinity-chromatography using immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human FES.
Concentration:	1 mg/ml
Gene Name:	FES
Other Name:	C-FES; FES/FPS; FPS; kinase Fes
Human Swiss Prot No.:	P07332
Formulation:	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Specificity:	FES Antibody antibody detects endogenous levels of FES.
Storage and Stability:	Stable at -20°C for at least 1 year.