

Antibodies Product Datasheet

Product Name:	STK39 Antibody
Catalog Number:	R12-3573
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	STK39 Antibody
Modification:	Unmodified/Total
Source:	Rabbit Polyclonal
Purification:	The antibody was purified from rabbit antiserum by affinity-chromatography using immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human STK39.
Concentration:	1 mg/ml
Gene Name:	STK39
Other Name:	STE20/SPS1-related proline-alanine-rich protein kinase; Ste-20-related kinase; Serine/threonine-protein kinase 39; DCHT; STK39; SPAK
Human Swiss Prot No.:	Q9UEW8
Formulation:	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Specificity:	STK39 Antibody antibody detects endogenous levels of STK39.
Storage and Stability:	Stable at -20°C for at least 1 year.