

Antibodies Product Datasheet

Product Name:	GPR120 Antibody
Catalog Number:	G127
Reactivity:	Human, Mouse, Rat
Applications:	WB, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human GPR120. AA range:221-270
Concentration:	1 mg/ml
Gene Name:	O3FAR1
Protein Name:	Omega-3 fatty acid receptor 1
Other Name:	O3FAR1; GPR120; GPR129; PGR4; Omega-3 fatty acid receptor 1; G-protein coupled receptor 120; G-protein coupled receptor 129; G-protein coupled receptor GT01; G-protein coupled receptor PGR4
Human Swiss Prot No.:	Q5NUL3
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	GPR120 Polyclonal Antibody detects endogenous levels of GPR120 protein.
Storage and Stability:	-20°C/1 year