

Antibodies Product Datasheet

Product Name:	PADI2 Antibody
Catalog Number:	C48016
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human PADI2
Concentration:	1 mg/ml
Gene Name:	PADI2 KIAA0994 PDI2
Protein Name:	Protein-arginine deiminase type-2 (EC 3.5.3.15) (PAD-H19) (Peptidylarginine deiminase II) (Protein-arginine deiminase type II)
Other Name:	
Human Swiss Prot No.:	Q9Y2J8
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of PADI2 at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year