

Antibodies Product Datasheet

Product Name:	Cleaved-Caspase-6 p11 (A194) Antibody
Catalog Number:	L0160
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000, ELISA 1:10000-20000
Modification:	Cleaved Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from Cleaved-Caspase-6 p11 (A194) . at AA range: 150-230
Concentration:	1 mg/ml
Gene Name:	CASP6 MCH2
Protein Name:	CleavedCaspase6p11
Other Name:	caspace 6, apoptosis-related cysteine peptidase
Human Swiss Prot No.:	P55212
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Cleaved-Caspase-6 p11 (A194) Polyclonal Antibody detects endogenous levels of Cleaved-Caspase-6 p11
Storage and Stability:	-20°C/1 year