

Antibodies Product Datasheet

Product Name:	ERCC2 Antibody
Catalog Number:	C43635
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 220-300
Concentration:	1 mg/ml
Gene Name:	ERCC2 XPD XPDC
Protein Name:	TFIIH basal transcription factor complex helicase XPD subunit (EC 3.6.4.12) (Basic transcription factor 2 80 kDa subunit) (BTF2 p80) (CXPB) (DNA excision repair protein ERCC-2) (DNA repair protein complementing XP-D cells) (TFIIH basal transcription facto
Other Name:	
Human Swiss Prot No.:	P18074
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	ERCC2 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year