

Antibodies Product Datasheet

Product Name:	SLP-76 Antibody
Catalog Number:	C16523
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. IHC-p: 1:100-300 ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from SLP-76 . at AA range: 270-350
Concentration:	1 mg/ml
Gene Name:	LCP2
Protein Name:	Lymphocyte cytosolic protein 2
Other Name:	LCP2; Lymphocyte cytosolic protein 2; SH2 domain-containing leukocyte protein of 76 kDa; SLP-76 tyrosine phosphoprotein; SLP76
Human Swiss Prot No.:	Q13094
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	SLP-76 Polyclonal Antibody detects endogenous levels of SLP-76 protein.
Storage and Stability:	-20°C/1 year