

Antibodies Product Datasheet

Product Name:	Anti-Human ZAP-70 (ABT246) Mouse Monoclonal Antibody
Catalog Number:	V0504
Clone ID:	ABT246
Reactivity:	Human, Mouse, Rat
Applications:	IHC, WB
Application Dilutions	IHC 1:100-200, WB 1:500-1000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
UniProt:	ZAP70
Gene Name:	P43403
Other Name:	70 kDa zeta associated protein, 70 kDa zeta-associated protein, EC 2.7.10.2, FLJ17670, FLJ17679, Selective T cell defect, SRK, STD, Syk related tyrosine kinase, Syk-related tyrosine kinase, Truncated ZAP kinase, Tyrosine protein kinase ZAP70, Tyrosine-protein kinase ZAP-70, TZK, ZAP 70, ZAP70, ZAP70_HUMAN, Zeta chain associated protein kinase 70kD, Zeta chain associated protein kinase 70 kDa, Zeta chain associated protein kinase 70 kDa isoform 1, Zeta chain associated protein kinase 70 kDa isoform 2, Zeta chain of T cell receptor associated protein kinase 70, Zeta chain TCR associated protein kinase 70kD, Zeta chain TCR associated protein kinase 70 kDa
Isotype	IgG1, kappa
Localization:	Tonsil/ Appendix
Positive Tissue Control:	Cytoplasmic, Nuclear
Immunogen	Synthetic peptide
Purification	Immunogen affinity purified
Formulation:	PBS, pH7.4, 50% glycerol, 0.03%Proclin300
Specificity:	The antibody can specifically recognize human ZAP-70 protein. In western blotting of Jurkat cell lysate, the antibody can label a 70 kDa band corresponding to ZAP-70.
Storage and Stability:	-20°C/1 year