

Antibodies Product Datasheet

Product Name:	RYR1 Antibody
Catalog Number:	C49799
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human RYR1
Concentration:	1 mg/ml
Gene Name:	RYR1 RYDR
Protein Name:	Ryanodine receptor 1 (RYR-1) (RyR1) (Skeletal muscle calcium release channel) (Skeletal muscle ryanodine receptor) (Skeletal muscle-type ryanodine receptor) (Type 1 ryanodine receptor)
Other Name:	
Human Swiss Prot No.:	P21817
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of RYR1 at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year