

Antibodies Product Datasheet

Product Name:	RBP4 mouse mAb
Catalog Number:	N1043
Reactivity:	Human
Applications:	ELISA
Dilutions:	ELISA 1:10000-20000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Purified recombinant human RBP4 protein fragments expressed in E.coli.
Concentration:	1 mg/ml
Gene Name:	rbp4
Protein Name:	
Other Name:	OTTHUMP00000020114;OTTHUMP00000020115;OTTHUMP00000020116;Plasma retinol binding protein 4;Plasma retinol-binding protein;Plasma retinol-binding protein(1-176);PRBP;PRO2222;RBP;RBP4;RDCCAS;RET4_HUMAN;Retinol binding protein 4;Retinol binding protein 4 inte
Human Swiss Prot No.:	P02753
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	
Storage and Stability:	-20°C/1 year