

Antibodies Product Datasheet

Product Name:	HSBP1 Antibody
Catalog Number:	C50843
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human HSBP1
Concentration:	1 mg/ml
Gene Name:	HSBP1 HSF1BP
Protein Name:	Heat shock factor-binding protein 1 (Nasopharyngeal carcinoma-associated antigen 13) (NPC-A-13)
Other Name:	
Human Swiss Prot No.:	O75506
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of HSBP1 at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year