

Antibodies Product Datasheet

Product Name:	53BP1 Antibody
Catalog Number:	B8454
Reactivity:	Human, Rat, Mouse
Applications:	IHC-p, IF(paraffin section), ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human TP53BP1. AA range:1744-1793
Concentration:	1 mg/ml
Gene Name:	TP53BP1
Protein Name:	Tumor suppressor p53-binding protein 1
Other Name:	TP53BP1; Tumor suppressor p53-binding protein 1; 53BP1; p53-binding protein 1; p53BP1
Human Swiss Prot No.:	Q12888
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	53BP1 Polyclonal Antibody detects endogenous levels of 53BP1 protein.
Storage and Stability:	-20°C/1 year