

Antibodies Product Datasheet

Product Name:	EIF4A2 Antibody
Catalog Number:	C48284
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human EIF4A2
Concentration:	1 mg/ml
Gene Name:	EIF4A2 DDX2B EIF4F
Protein Name:	Eukaryotic initiation factor 4A-II (eIF-4A-II) (eIF4A-II) (EC 3.6.4.13) (ATP-dependent RNA helicase eIF4A-2)
Other Name:	
Human Swiss Prot No.:	Q14240
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of EIF4A2 at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year