

Antibodies Product Datasheet

Product Name:	FoxL2 (Phospho-Ser263) Antibody
Catalog Number:	A1202
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human FOXL2 around the phosphorylation site of Ser263. AA range:229-278
Concentration:	1 mg/ml
Gene Name:	FOXL2
Protein Name:	Forkhead box protein L2
Other Name:	FOXL2; Forkhead box protein L2
Human Swiss Prot No.:	P58012
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Phospho-FoxL2 (S263) Polyclonal Antibody detects endogenous levels of FoxL2 protein only when phosphorylated at S263.
Storage and Stability:	-20°C/1 year