

Antibodies Product Datasheet

Product Name:	Cbl (Phospho-Tyr731) Antibody
Catalog Number:	P12-1005
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Cbl (Phospho-Tyr731) Antibody
Modification:	Phospho Specific
Source:	Rabbit Polyclonal
Purification:	The antibody was purified from rabbit antiserum by affinity-chromatography using phospho peptide. The antibody against non-phospho peptide was removed by chromatography using corresponding non-phospho peptide.
Immunogen:	The antiserum was produced against synthesized phospho-peptide derived from human Cbl around the phosphorylation site of Tyr731.
Concentration:	1 mg/ml
Gene Name:	CBL
Other Name:	E3 ubiquitin-protein ligase CBL; Casitas B-lineage lymphoma proto-oncogene; Proto-oncogene c-Cbl; RING finger protein 55; Signal transduction protein CBL
Human Swiss Prot No.:	P22681
Formulation:	Rabbit IgG in phosphate buffered saline (Phospho-without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Specificity:	Cbl (Phospho-Tyr731) antibody detects endogenous levels of Cbl only when phosphorylated at Tyr731.
Storage and Stability:	Stable at -20°C for at least 1 year.