

Antibodies Product Datasheet

Product Name:	NFAT5 Antibody
Catalog Number:	B1220
Reactivity:	Human, Mouse
Applications:	IHC-p, IF(paraffin section), ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human NFAT5. AA range:1171-1220
Concentration:	1 mg/ml
Gene Name:	NFAT5
Protein Name:	Nuclear factor of activated T-cells 5
Other Name:	NFAT5; KIAA0827; TONEBP; Nuclear factor of activated T-cells 5; NFAT5; T-cell transcription factor NFAT5; Tonicity-responsive enhancer-binding protein; TonE-binding protein; TonEBP
Human Swiss Prot No.:	O94916
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	NFAT5 Polyclonal Antibody detects endogenous levels of NFAT5 protein.
Storage and Stability:	-20°C/1 year