

## Antibodies Product Datasheet

<b>Product Name:</b>	AI-BP Antibody
<b>Catalog Number:</b>	C14522
<b>Reactivity:</b>	Human, Mouse, Rat
<b>Applications:</b>	WB, IF/ICC, ELISA
<b>Dilutions:</b>	Western Blot: 1/500 - 1/2000. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/40000. Not yet tested in other applications.
<b>Modification:</b>	Unmodified/Total
<b>Source:</b>	Polyclonal Rabbit
<b>Purification:</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Immunogen:</b>	The antiserum was produced against synthesized peptide derived from human AIBP. AA range:121-170
<b>Concentration:</b>	1 mg/ml
<b>Gene Name:</b>	APOA1BP
<b>Protein Name:</b>	NAD(P)H-hydrate epimerase
<b>Other Name:</b>	APOA1BP; AIBP; YJEFN1; NAD(P)H-hydrate epimerase; Apolipoprotein A-I-binding protein; AI-BP; NAD(P)HX epimerase; YjeF N-terminal domain-containing protein 1; YjeF_N1
<b>Human Swiss Prot No.:</b>	Q8NCW5
<b>Formulation:</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Specificity:</b>	AI-BP Polyclonal Antibody detects endogenous levels of AI-BP protein.
<b>Storage and Stability:</b>	-20°C/1 year