

Antibodies Product Datasheet

Product Name:	Cdc50B Antibody
Catalog Number:	C12119
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human TMEM30B. AA range:91-140
Concentration:	1 mg/ml
Gene Name:	TMEM30B
Protein Name:	Cell cycle control protein 50B
Other Name:	TMEM30B; CDC50B; Cell cycle control protein 50B; Transmembrane protein 30B
Human Swiss Prot No.:	Q3MIR4
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Cdc50B Polyclonal Antibody detects endogenous levels of Cdc50B protein.
Storage and Stability:	-20°C/1 year